

THE NEOTIA UNIVERSITY

WORK INSTRUCTION

1.0 EXPERIMENT NO: BO/01

2.0 NAME OF EXPERIMENT: DISPERSIVE POWER

3.0 OBJECTIVE: DETERMINATION OF DISPERSIVE POWER OF THE MATERIAL OF A PRISM

4.0 **PRINCIPLE:** If A is the angle of the prism and δ_m is the minimum deviation of a ray of monochromatic light refracted through the prism in a principle section, then the refractive index of the material of the prism μ_λ for the light of the given wavelength is expressed by

$$\mu_\lambda = \frac{\sin \frac{(A + \delta_m)}{2}}{\sin \frac{A}{2}} \dots\dots\dots (1)$$


if μ_1 and μ_2 are the refractive indices of the material of the prism for lights of two given colours and μ is the refractive index of the material for the mean colour between the two given colours, then the dispersive power of the material of the prism ω between the wavelength region considered is given by

$$\omega = \frac{\mu_1 - \mu_2}{\mu - 1} \dots\dots\dots (2)$$

Equations (1) & (2) are used as the working formula for the determination of the dispersive power of the material of a prism.

5.0 TOOLS/APPARATUS REQUIRED:

- 5.1 Spectrometer
- 5.2 Mercury lamp
- 5.3 Prism
- 5.4 Sprit level



DISPERSIVE POWER SET - UP

THE NEOTIA UNIVERSITY

Procedure :

Using spirit level, level the base of the spectrometer, prism table, telescope and collimator.

↓
At the minimum deviation position of prism for yellow line, adjust the telescope and collimator for parallel rays

↓
Determine the vernier constant (V.C.) of the spectrometer

↓
Place the prism perpendicular to the collimator axis

↓
There will be two refracted images in two sides of the prism

↓
Touching the cross-wire on the image in both sides and take the reading of two vernier position
Determine the angle (A) of the prism from tab. 2

↓
Place the prism at the minimum deviation position and take the reading of two vernier position for three different colours

↓
Take direct reading of the telescope

↓
Determine the angle of minimum deviation from tab. – 5

↓
Determine the refractive index of three different colours

↓
Determine the dispersive power of the material of the prism using equn. – 2

Tabulation:

TABLE-1

Determination of vernier constant (V.C.) of the spectrometer

..... divisions(say m) of the vernier scale=..... divisions (say n) of the main scale

| Value of 1 smallest circular scale division (l_1) (min/sec) | Value of 1 division of the vernier scale $\frac{n}{m}$ (min/sec) | Vernier constant of the spectrometer v.c. $(1 - \frac{n}{m}) \times l_1$ (min/sec) |
|---|--|---|
| | | |

THE NEOTIA UNIVERSITY

TABLE-2
Determination of the angle (A) of the prism

| Vernier scale | Reading of the 1 st position of the telescope | | | Reading of the 2 nd position of the telescope | | | 2A = (a ~ b) (deg) | Mean 2A (deg) | Angle (A) of the prism (deg) |
|-------------------------|--|-------------------|---|--|-------------------|---|--------------------|---------------|------------------------------|
| | Main scale reading m_s (deg) | Vernier no. v_s | Total $m_s + v_s \times v.c$ (deg). (a) | Main scale reading m_s (deg) | Vernier no. v_s | Total $m_s + v_s \times v.c$ (deg). (b) | | | |
| 1 st Vernier | | | | | | | | | |
| 2 nd Vernier | | | | | | | | | |

TABLE-3
Reading of the telescope at the position of minimum deviation :

| Colour of the line | Vernier scale | No. of obs. | Main scale reading m_s (deg) | Vernier no. v_s | Total $m_s + v_s \times v.c$ (deg). | Mean reading (deg). |
|--------------------|-------------------------|-------------|--------------------------------|-------------------|-------------------------------------|---------------------|
| Red | 1 st Vernier | 1 | | | | d ₁ |
| | | 2 | | | | |
| | | 3 | | | | |
| | 2 nd Vernier | 1 | | | | d ₂ |
| | | 2 | | | | |
| | | 3 | | | | |
| Yellow | 1 st Vernier | 1 | | | | d ₃ |
| | | 2 | | | | |
| | | 3 | | | | |
| | 2 nd Vernier | 1 | | | | d ₄ |
| | | 2 | | | | |
| | | 3 | | | | |
| Blue | 1 st Vernier | 1 | | | | d ₅ |
| | | 2 | | | | |
| | | 3 | | | | |
| | 2 nd Vernier | 1 | | | | d ₆ |
| | | 2 | | | | |
| | | 3 | | | | |

THE NEOTIA UNIVERSITY

TABLE-4
Direct reading of the telescope:

| Vernier scale | No. of obs. | Main scale reading m_s (deg) | Vernier no. v_s | Total $m_s + v_s \times v.c$ (deg). | Mean reading (deg). |
|-------------------------|-------------|--------------------------------|-------------------|-------------------------------------|---------------------|
| 1 st Vernier | 1 | | | | d ₇ |
| | 2 | | | | |
| | 3 | | | | |
| 2 nd Vernier | 1 | | | | d ₈ |
| | 2 | | | | |
| | 3 | | | | |

TABLE-5
Determination of angle of minimum deviation:

| Colour of the line | Angle of minimum deviation δ_m (deg). | | Mean δ_m (deg). |
|--------------------|--|---------------------------------|------------------------|
| | 1 st Vernier | 2 nd Vernier | |
| Red | d ₁ ~ d ₇ | d ₂ ~ d ₈ | |
| Yellow | d ₃ ~ d ₇ | d ₄ ~ d ₈ | |
| Blue | d ₅ ~ d ₇ | d ₆ ~ d ₈ | |

TABLE-6
Determination of dispersive power of the material of the prism:

| Angle (A) of the prism (deg) | Colour of the line | δ_m (deg). | Refractive index | Dispersive power $\omega = (\mu_1 - \mu_2) / (\mu - 1)$ |
|------------------------------|--------------------|-------------------|------------------|---|
| | Red | | μ_1 | |
| | Yellow | | μ | |
| | Blue | | μ_2 | |

Discussion:

You have to write all the difficulties you faced during the experiment and their remedies. Also you have to mention some way out that one should adopt during the practical to have a better result.

Viva voice: go through the chapter of young's modulus and elasticity from these books.

- 1) OPTICS - Ghatak
- 2) OPTICS – K. G. Majumdar
- 3) ADVANCED PRACTICAL PHYSICS- Ghosh & Majumdar
- 4) PRACTICAL PHYSICS- Rakshit, Chatterjee & Saha